

# SARS Nucleocapsid Recombinant Protein



RPPB5318

---

## Product Information    Protein Information

**Product SKU:**

RPPB5318

**Accession:**

P0C6X7.1

**Host:**

Escherichia Coli.

**Protein description:**

The E.coli derived recombinant protein contains the Nucleocapsid protein 1-49 amino acids immunodominant regions.

**Appearance:**

Sterile Filtered clear solution.

**Formulation:**

SARS Nucleocapsid protein solution is supplied in PBS.

**Purity:**

Protein is >90% pure as determined SDS-PAGE.

**Stability:**

Protein is shipped on ice packs. Upon arrival, Store at -20°C.