

SARS Nucleocapsid (340-390) Recombinant Protein



RPPB5319

Product Information Protein Information

Product SKU:

RPPB5319

Accession:

POC6X7.1

Host:

Escherichia Coli.

Protein description:

The E.coli derived recombinant protein contains the Nucleocapsid C-Terminal protein 340-390 amino acids immunodominant regions fused to 6xHis tag at C-terminal.

Appearance:

Sterile Filtered clear solution.

Formulation:

SARS Nucleocapsid protein solution is supplied in PBS.

Purity:

Protein is >90% pure as determined SDS-PAGE.

Stability:

Protein is shipped on ice packs. Upon arrival, Store at -20°C.