

SAR-CoV-2 CoV-2-S1 (319-541) Recombinant Protein



RPPB5330

Product Information Protein Information

Product SKU:

RPPB5330

Accession:

QHD43416.1

Host:

HEK293 Cells.

Protein description:

The HEK293 derived recombinant protein contains the Coronavirus 2019 Spike Glycoprotein S1 Receptor Binding Domain [RBD], Wuhan-Hu-1 strain, amino acids 319-541 fused to Fc tag at C-terminal.

Appearance:

Lyophilized freeze-dried powder.

Formulation:

CoV-S1 RBD protein is lyophilized from 1x PBS pH-7.4 + 5% trehalose.

Purity:

Protein is >95% pure as determined SDS-PAGE.

Stability:

Lyophilized Cov-2 S1 Glycoprotein RBD although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution CoV2 Spike protein should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.