

SAR-CoV-2 CoV-2 S1 (16-685) Recombinant Protein



RPPB5335

Product Information Protein Information

Product SKU:

RPPB5335

Accession:

QHD43416.1

Host:

HEK293 Cells.

Protein description:

The HEK293 derived Biotinylated recombinant protein contains the Coronavirus 2019 CoV-2 Spike Glycoprotein S1, Wuhan-Hu-1 strain, amino acids 16-685 fused to His tag & Avi-tag at C-terminal.

Appearance:

Sterile Filtered White lyophilized (freeze-dried) powder.

Formulation:

CoV-2 S1 protein is supplied in 1x PBS pH-7.4 & 5% trehalose.

Purity:

Protein is >90% pure as determined SDS-PAGE.

Stability:

Lyophilized Cov-2 Spike S1 Glycoprotein although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution CoV2 Spike protein should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.