

SAR-CoV-2 CoV-2 S2 Recombinant Protein



RPPB5336

Product Information Protein Information

Product SKU:

RPPB5336

Accession:

YP_009724390.1

Host:

HEK293

Protein description:

The HEK293 derived recombinant protein contains the Novel Coronavirus 2019-nCoV Spike Glycoprotein S2, Wuhan-Hu-1 strain, amino acids 685-1211 fused to Fc tag at C-terminal.

Appearance:

Sterile Filtered clear solution.

Formulation:

nCoV-S2 protein solution is supplied in DPBS.

Purity:

Protein is >85% pure as determined SDS-PAGE.

Stability:

Protein is shipped on ice packs. Upon arrival, Store at -20°C.