

SAR-CoV-2 CoV-2 S2 Recombinant Protein



RPPB5337

Product Information Protein Information

Product SKU:

RPPB5337

Accession:

YP_009724390.1

Host:

Sf9, Baculovirus Cells.

Protein description:

The Sf9 derived recombinant protein contains the Coronavirus 2019 CoV-2 Spike Glycoprotein S2, Wuhan-Hu-1 strain, 685-1211 amino acids, having a Mw of 60.1 kDa and fused to 6xHis tag at C-terminal.

Appearance:

Sterile Filtered clear solution.

Formulation:

CoV-2 S2 protein solution is supplied in DPBS.

Purity:

Protein is >85% pure as determined SDS-PAGE.

Stability:

CoV-2 S2 Glycoprotein is shipped on ice packs. Upon arrival, Store at -20°C.