

SAR-CoV-2 CoV-2 Spike (300-600) Recombinant Protein



RPPB5341

Product Information Protein Information

Product SKU:

RPPB5341

Accession:

NC_045512.2

Host:

Escherichia Coli.

Protein description:

The E.Coli derived recombinant protein contains the Coronavirus 2019 Spike Receptor Binding Domain (300-600 a.a.) immunodominant region, fused to 6xHis tag at C-terminal

Appearance:

Sterile Filtered clear solution.

Formulation:

CoV 2019 Spike Receptor Binding Domain Protein 1mg/ml solution is supplied in 1x PBS.

Purity:

CoV 2019 Spike Receptor Binding Domain Protein is >90% pure as determined SDS-PAGE.

Stability:

CoV 2019 Spike Receptor Binding Domain Protein is shipped on ice packs. Upon arrival, Store at -20°C.