

SAR-CoV-2 CoV-2-Spike (1-260) Recombinant Protein



RPPB5342

Product Information Protein Information

Product SKU:

RPPB5342

Accession:

YP_009724390.1

Host:

HEK293

Protein description:

The HEK293 derived recombinant protein contains the Coronavirus 2019-Spike N-Terminal Domain, Wuhan-Hu-1 strain, amino acids 1-260 fused to Fc tag at C-terminal.

Appearance:

Sterile Filtered clear solution.

Formulation:

CoV-2 spike protein solution is supplied in DPBS.

Purity:

Protein is >85% pure as determined SDS-PAGE.

Stability:

CoV-2 Spike Protein is shipped on ice packs. Upon arrival, Store at -20°C.