

SAR-CoV-2 CoV-2 Membrane Env. Recombinant Protein



RPPB5344

Product Information Protein Information

Product SKU:

RPPB5344

Host:

E.Coli

Protein description:

The E.Coli derived recombinant fusion protein contains the Coronavirus 2019 Full-Length Membrane and Envelope proteins, Wuhan-Hu-1 strain, having a Mw of 34.2 kDa fused to His tag at C-terminal.

Appearance:

Sterile Filtered clear solution.

Formulation:

CoV-2 Membrane Envelope protein solution is supplied in 1x PBS.

Purity:

Protein is >90% pure as determined SDS-PAGE.

Stability:

CoV-2 Membrane Envelope fusion Protein is shipped on ice packs. Upon arrival, Store at -20°C.