

SAR-CoV-2 CoV-2 N (201-419) Recombinant Protein



RPPB5351

Product Information Protein Information

Product SKU:

RPPB5351

Accession:

QHO62115.1

Host:

Escherichia Coli.

Protein description:

The E.Coli derived recombinant protein contains the Coronavirus 2019 C-Terminal domain nucleopocapsid (201-419 a,a.) fused to 6xHis tag at C-terminal and having a molecular weight of 23.9kDa.

Appearance:

Sterile Filtered clear solution.

Formulation:

CoV 2019 Nucleocapsid-c-terminal domain Protein solution is supplied in 1x PBS.

Purity:

CoV 2019 Nucleocapsid-c-terminal Protein ein is >90% pure as determined SDS-PAGE.

Stability:

CoV 2019 Nucleocapsid c-terminal domain Protein is shipped on ice packs. Upon arrival, Store at -20°C.