

Borrelia p41 Recombinant Protein



RPPB5352

Product Information Protein Information

Product SKU:

RPPB5352

Host:

Escherichia Coli.

Protein description:

The E.Coli derived recombinant protein contains the p41 immunodominant regions, 158-296 amino acids. Borrelia p41 is fused to 6xHis tag at N-terminal and purified by proprietary chromatographic techniques.

Formulation:

25mM glycine, pH 9.6 and 50% glycerol.

Purity:

Protein is >90% pure as determined by SDS-PAGE.

Stability:

Borrelia p41 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.