

Borrelia Miyamotoi GlpQ Recombinant Protein



RPPB5376

Product Information Protein Information

Product SKU:

RPPB5376

Host:

Escherichia Coli.

Protein description:

Recombinant Borrelia Miyamotoi GlpQ produced in E.coli is a non-glycosylated, polypeptide chain having a calculated molecular mass of 38,132 Dalton. Borrelia Miyamotoi GlpQ is expressed with a 10xHis tag at N-terminus and purified by proprietary chromatographic techniques.

Appearance:

Sterile Filtered clear solution.

Formulation:

Borrelia Miyamotoi GlpQ is supplied in 20mM HEPES buffer pH-7.6, 250mM NaCl and 20% glycerol.

Purity:

Greater than 80.0% as determined by SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.