

# HSV-1 gD Recombinant Protein



RPPB5384

---

## Product Information    Protein Information

**Product SKU:**

RPPB5384

**Protein description:**

The E.Coli derived recombinant protein contains the HSV-1 gD immunodominant regions, 266-394 amino acids and fused to a GST-Tag at C-terminus.

**Appearance:**

Sterile Filtered clear solution.

**Formulation:**

25mM Tris-HCl pH 8, 1mM EDTA, and 50% glycerol.

**Purity:**

HSV-1 gD protein is >95% pure as determined by 10% PAGE (coomassie staining).

**Stability:**

HSV-1 gD protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.