

HSV-1 gD (84-175) Recombinant Protein



RPPB5385

Product Information Protein Information

Product SKU:

RPPB5385

Protein description:

The E.Coli derived recombinant protein contains the HSV-1 gD (84-175) immunodominant region.

Host:

Escherichia Coli.

Appearance:

Sterile Filtered clear solution.

Formulation:

0.1 % SDS, 50% glycerol and 100 mM NaCl.

Purity:

Protein is >95% pure as determined by SDS-PAGE

Stability:

HSV-1 gD although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.