

# HSV-1 gG Recombinant Protein



RPPB5386

---

## Product Information    Protein Information

**Product SKU:**

RPPB5386

**Protein description:**

The E.Coli derived recombinant protein contains the HSV-1 gG immunodominant regions, 84-175 amino acids and fused to a GST-Tag at C-terminus.

**Appearance:**

Sterile Filtered clear solution.

**Formulation:**

25mM Tris-HCl pH 7.2, 1mM EDTA, and 50% glycerol.

**Purity:**

HSV-1 gG protein is >95% pure as determined by 10% PAGE (coomassie staining).

**Stability:**

HSV-1 gG protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.