

HSV-2 gD Recombinant Protein



RPPB5387

Product Information Protein Information

Product SKU:

RPPB5387

Protein description:

The E.Coli derived 39.7 kDa recombinant protein contains the HSV-2 gD immunodominant regions 266-394 amino acids, fused with 26 kDa GST-tag.

Appearance:

Sterile Filtered clear solution.

Formulation:

25mM Tris pH 7.2, 1mM EDTA, and 50% glycerol.

Purity:

HSV-2 gD protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

HSV-2 gD Protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.