

HSV-2 gD (525-578) Recombinant Protein



RPPB5388

Product Information Protein Information

Product SKU:

RPPB5388

Host:

Escherichia Coli.

Protein description:

The E.Coli derived recombinant protein contains the HSV-2 gD (525-578) immunodominant region. HSV2 gD is fused to a six histidine his tag at c-terminus and purified by proprietary chromatographic techniques.

Appearance:

Sterile Filtered clear solution.

Formulation:

(1mg/ml) in 1x PBS

Purity:

Protein is >95% pure as determined by SDS-PAGE

Stability:

HSV-2 gD although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.