

HSV-2 gG Recombinant Protein



RPPB5390

Product Information Protein Information

Product SKU:

RPPB5390

Protein description:

The E.Coli derived recombinant protein contains the HSV-2 gG immunodominant regions, 525-578 amino acids, the total MW is 32,100 Dalton, fused with 26 kDa GST-tag.

Appearance:

Sterile Filtered clear solution.

Formulation:

25mM Tris pH 8, 1.5M Urea and 50% glycerol.

Purity:

HSV-2 gG protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

HSV-2 gG protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.