

HSV-8 M Recombinant Protein



RPPB5391

Product Information Protein Information

Product SKU:

RPPB5391

Protein description:

The E.Coli derived recombinant protein contains the C-terminal immunodominant regions from ORF65 140-170 a.a. and N-terminal regions from ORF8 32-62 a.a. The protein is fused with a GST tag.

Appearance:

Sterile Filtered clear solution.

Formulation:

100mM NaCl, 0.1% SDS and 50% glycerol.

Purity:

HSV-8 Mosaic protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

HSV-8 Mosaic protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.