

Chlamydia W2 Recombinant Protein



RPPB5395

Product Information Protein Information

Product SKU:

RPPB5395

Host:

Escherichia Coli.

Protein description:

The E.coli derived recombinant 6x His fusion at C-terminus protein contains Chlamydia Trachomatis MOMP protein epitopes, 66-165 amino acids.

Formulation:

20mM Tris-HCl, pH 7.2, 50% glycerol and 1.5M urea.

Purity:

Protein is >90% pure as determined by SDS-PAGE.

Stability:

Chlamydia W2 protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.