

Chlamydia W4 Recombinant Protein



RPPB5396

Product Information Protein Information

Product SKU:

RPPB5396

Host:

Escherichia Coli.

Protein description:

The E.coli derived recombinant 6xHis fusion at C-terminus protein contains Chlamydia Trachomatis MOMP protein epitopes, 191-286 amino acids.

Formulation:

(1mg/1ml) 20mM Tris-HCl, pH 7.2, 50% glycerol and 1.5M urea.

Purity:

Protein is >90% pure as determined by SDS PAGE.

Stability:

Chlamydia W4 protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.