

Chlamydia W5 Recombinant Protein



RPPB5397

Product Information Protein Information

Product SKU:

RPPB5397

Protein description:

The E.coli derived recombinant 6x His fusion at C-terminus protein contains Chlamydia Trachomatis MOMP protein epitopes, 252-354 amino acids, and a 6xHis Tag fused at C-terminus.

Formulation:

20mM Tris-HCl, pH 7.2, 1.5M urea and 50% glycerol.

Purity:

Protein is >90% pure as determined by SDS- PAGE.

Stability:

Chlamydia W5 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.