

# Chlamydia W4-W5 Recombinant Protein



RPPB5398

---

## Product Information    Protein Information

**Product SKU:**

RPPB5398

**Host:**

Escherichia Coli.

**Protein description:**

The E.coli derived recombinant 6xHis fusion at C-terminus protein contains Chlamydia Trachomatis MOMP protein epitopes, 191-354 amino acids.

**Formulation:**

10mM Tris-HCl, pH 6.5, 100mM Sodium Phosphate and 8M urea.

**Purity:**

Chlamydia W4-W5 protein is >95% pure as determined by 10% PAGE (coomassie staining) and RP-HPLC.

**Stability:**

Chlamydia W4-W5 protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.