

# Chlamydia W5-W6 Recombinant Protein



RPPB5399

---

## Product Information    Protein Information

**Product SKU:**

RPPB5399

**Host:**

Escherichia Coli.

**Protein description:**

The E.coli derived recombinant 6xHis fusion at C-terminus protein contains Chlamydia Trachomatis MOMP protein epitopes, 252-398 amino acids.

**Formulation:**

20mM Tris-HCl, pH 7.2, 50% glycerol and 1.5M urea (1mg/1ml).

**Purity:**

Protein is >90% pure as determined by SDS-PAGE.

**Stability:**

Chlamydia W5-W6 protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.