

Chlamydia W3-W6 Recombinant Protein



RPPB5400

Product Information Protein Information

Product SKU:

RPPB5400

Protein description:

The E.coli derived recombinant 6xHis fusion at C-terminus. protein contains Chlamydia TrachomatisMOMP protein epitopes, 128-398 amino acids.

Formulation:

10mM Tris-HCl, pH 6.5, 100mM Sodium Phosphate and 8M urea.

Purity:

Chlamydia W3-W6 protein is >95% pure as determined by 10% PAGE (coomassie staining) and RP-HPLC.

Stability:

Chlamydia W3-W6 protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.