

# Chlamydia PGP-3D Recombinant Protein



RPPB5401

---

## Product Information   Protein Information

**Product SKU:**

RPPB5401

**Host:**

Escherichia Coli.

**Protein description:**

The E.coli derived recombinant protein contains Chlamydia Trachomatis PGP-3D full length protein epitope. The protein is fused to a 6xHis-Tag at C-terminus.

**Formulation:**

10mM Tris-HCl, pH 8, 1.5M Urea and 50% Glycerol.

**Purity:**

Protein is >90% pure as determined by SDS-PAGE.

**Stability:**

Chlamydia PGP-3D although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.