

# Chlamydia HSP70 (462-503 a.a.) Recombinant Protein



RPPB5402

---

## Product Information    Protein Information

**Product SKU:**

RPPB5402

**Host:**

Escherichia Coli.

**Protein description:**

The E.coli derived recombinant protein contains Chlamydia TrachomatisHSP70 protein epitopes, 462-503 amino acids.

**Formulation:**

50mM Tris-HCl, pH 8.0, 8M Urea, 60 mM NaCl and 50% glycerol.

**Purity:**

Chlamydia HSP70 protein is >95% pure as determined by 10% PAGE (coomassie staining) and RP-HPLC.

**Stability:**

Chlamydia HSP70 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.