

Chlamydia HSP70 (549-660a.a.) Recombinant Protein



RPPB5403

Product Information Protein Information

Product SKU:

RPPB5403

Host:

Escherichia Coli.

Protein description:

The E.coli derived recombinant protein contains Chlamydia TrachomatisHSP70 protein epitopes, 549-660 amino acids.

Formulation:

50mM Tris-HCl, pH 8.0, 8M Urea, 60mM NaCl and 50% glycerol.

Purity:

Chlamydia HSP70 protein is >95% pure as determined by 10% PAGE (coomassie staining) and RP-HPLC.

Stability:

Chlamydia HSP70 protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.