

HCV Core Genotype-1a Recombinant Protein



RPPB5406

Product Information Protein Information

Product SKU:

RPPB5406

Protein description:

The E.coli derived recombinant protein fused to a His tag at C-terminus contains the HCV core nucleocapsid immunodominant regions, amino acids 2-119.

Formulation:

1.5M urea, 25mM Tris-HCl pH 8.0, 0.2% Triton-X and 50% Glycerol.

Purity:

HCV Core Genotype-1a protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

HCV Core Genotype-1a although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.