

# HCV Core Genotype-3/10 Recombinant Protein



RPPB5416

---

## Product Information    Protein Information

**Product SKU:**

RPPB5416

**Host:**

Escherichia Coli.

**Protein description:**

The E.coli derived recombinant protein contains the HCV core nucleocapsid immunodominant regions, amino acids 2-119.

**Formulation:**

50mM Tris-HCl, pH-8, 60mM NaCl, 10mM glutathione, 0.25% sarkosil & 50% glycerol.

**Purity:**

HCV Core Genotype-3/10 protein is >95% pure as determined by 10% PAGE (coomassie staining).

**Stability:**

HCV Core Genotype-3/10 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.