

# HCV Core Genotype-1b Biotin Recombinant Protein



RPPB5418

---

## Product Information    Protein Information

**Product SKU:**

RPPB5418

**Protein description:**

The E.coli derived recombinant Biotin Labeled protein contains the HCV core nucleocapsid immunodominant regions, amino acids 2-119, having an MW of 22kDa. The protein is fused to a beta-galactosidase (114 kDa) at the N-terminus.

**Formulation:**

20mM Tris-HCl pH 8 and 8M urea.

**Purity:**

Protein is >95% pure as determined by SDS-PAGE.

**Stability:**

HCV Core although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.