

# HCV Core Recombinant Protein



RPPB5419

---

## Product Information    Protein Information

**Product SKU:**

RPPB5419

**Protein description:**

The E.coli derived recombinant HRP Labeled protein contains the HCV core nucleocapsid immunodominant regions, amino acids 2-119.

**Formulation:**

25mM Tris-HCl pH-8, 1mM EDTA, 1.5M urea and 50%glycerol.

**Purity:**

HCV Core protein is >95% pure as determined by 10% PAGE (coomassie staining).

**Stability:**

HCV Core HRP although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.