

# HCV Core NS3 Recombinant Protein



RPPB5421

---

## Product Information    Protein Information

**Product SKU:**

RPPB5421

**Host:**

Escherichia Coli.

**Protein description:**

Recombinant HCV protein produced in E.Coli containing the Core (130 a.a.) and NS3 (236 a.a.) excluding NS4 and NS5 domains having a total Mw of 45kDa and purified by proprietary chromatographic techniques.

**Appearance:**

Sterile filtered colorless solution.

**Formulation:**

Sterile filtered solution containing PBS, 1M urea and 25mM Arginine.

**Purity:**

Greater than 95.0% as determined by SDS-PAGE.

**Stability:**

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.