

HCV Core 22kDa Recombinant Protein



RPPB5427

Product Information Protein Information

Product SKU:

RPPB5427

Protein description:

The E.coli derived recombinant protein contains the HCV core nucleocapsid immunodominant regions, amino acids 2-192. The FLSCN labeled protein is fused with b-galactosidase (114 kDa) at N-terminus.

Formulation:

20mM Tris-HCl pH 8, 8M urea and 10mM B-ME.

Purity:

HCV-Core protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

HCV-Core although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.