

# HCV Core 22kDa Recombinant Protein



RPPB5428

---

## Product Information    Protein Information

**Product SKU:**

RPPB5428

**Protein description:**

The E.coli derived recombinant protein contains the HCV core nucleocapsid immunodominant regions, amino acids 2-192. The Rhodamine labeled protein is fused with b-galactosidase (114 kDa) at N-terminus.

**Formulation:**

20mM Tris-Hcl pH 8, 8M urea and 10mM B-ME.

**Purity:**

HCV-Core protein is >95% pure as determined by 10% PAGE (coomassie staining).

**Stability:**

HCV-Core although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.