

HCV NS3 1a Recombinant Protein



RPPB5433

Product Information Protein Information

Product SKU:

RPPB5433

Protein description:

The E.coli derived recombinant protein contains the HCV NS3 immunodominant regions, amino acids 1356-1459. The protein is fused to a GST tag at N-Terminus.

Formulation:

1.5M urea, 25mM Tris-HCl pH-8, 0.2% Triton-X & 50% Glycerol.

Purity:

HCV NS3 Genotype-1a protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

HCV NS3 Genotype-1a although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.