

# HCV NS3 (1192-1459 a.a.) Recombinant Protein



RPPB5434

---

## Product Information    Protein Information

**Product SKU:**

RPPB5434

**Protein description:**

The E.coli derived recombinant protein contains the HCV NS3 immunodominant regions, amino acids 1192-1459. The protein is fused to a GST tag at N-Terminus.

**Formulation:**

1.5M urea, 25mM Tris-HCl pH-8, 0.2% Triton-X & 50% Glycerol.

**Purity:**

HCV NS3 Genotype-1b protein is >95% pure as determined by 10% PAGE (coomassie staining).

**Stability:**

NS3 Genotype-1b although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.