

HCV NS3 Genotype-1b Recombinant Protein



RPPB5435

Product Information Protein Information

Product SKU:

RPPB5435

Protein description:

The E.coli derived 26.2 kDa recombinant protein contains the HCV NS3 c33c immunodominant regions. The 252 amino acid protein is fused to 6xHis tag.

Formulation:

50mM NaPO₄ pH 8.3 and 10mM DTT.

Purity:

HCV NS3 Genotype-1b protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

HCV NS3 Genotype-1b although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.

Amino Acid Sequence:

Starts: MRDSQSQTFFQ Ends: DSVDCNTCVT.