

# HCV NS3 Genotype-1c Recombinant Protein



RPPB5437

---

## Product Information    Protein Information

**Product SKU:**

RPPB5437

**Protein description:**

The E.coli derived recombinant His-Tag fusion protein contains the HCV NS3 immunodominant regions, amino acids 1192-1459.

**Formulation:**

1.5M urea, 25mM Tris-HCl pH-8, 0.2% Triton-X & 50% Glycerol.

**Purity:**

HCV NS3 Genotype-1c protein is >95% pure as determined by 10% PAGE (coomassie staining).

**Stability:**

HCV NS3 Genotype-1c although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.