

HCV NS3 Genotype-3 Recombinant Protein



RPPB5442

Product Information Protein Information

Product SKU:

RPPB5442

Protein description:

The E.coli derived recombinant protein contains the HCV NS3 immunodominant regions, amino acids 1359-1456. The protein is fused to a GST tag at N-Terminus.

Formulation:

1.5M urea, 25mM Tris-HCl pH-8, 0.2% Triton-X & 50% Glycerol.

Purity:

HCV NS3 Genotype-3 protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

HCV NS3 Genotype-3 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.