

HCV NS3 Genotype-5a Recombinant Protein



RPPB5446

Product Information Protein Information

Product SKU:

RPPB5446

Protein description:

E.coli derived recombinant from HCV subtype 5a fused to a GST tag. The protein contains the full-length HCV NS3 (c33c) immunodominant region.

Formulation:

25mM Tris-HCl, 1mM EDTA, 1.5M urea & 50% glycerol.

Purity:

HCV NS3 Genotype-5a protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

HCV NS3 Genotype-5a although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.