

# HCV NS3 Recombinant Protein



RPPB5449

---

## Product Information    Protein Information

**Product SKU:**

RPPB5449

**Protein description:**

The E.coli derived 22kDa recombinant protein contains the HCV NS3 genotype 1b immunodominant regions, amino acids 1450-1643. The protein is fused with b-galactosidase (114 kDa) at N-terminus, pI 5.43.

**Formulation:**

20mM Tris-HCl pH-8, 8M urea & 10mM B-ME.

**Purity:**

HCV NS3 protein is >95% pure as determined by 10% PAGE (coomassie staining).

**Stability:**

HCV NS3 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.