

HCV NS3 Recombinant Protein



RPPB5451

Product Information Protein Information

Product SKU:

RPPB5451

Protein description:

The E.coli derived recombinant Biotin Labeled protein contains the HCV NS3 immunodominant regions, a.a. 1450-1643, 22 kDa. The Biotin labeled protein is fused with a 6xHis-Tag at N-terminus.

Formulation:

(1mg/ml) 1.5M urea, 20mM Tris-HCl pH 8.0 and 10mM β -mercaptoethanol.

Purity:

HCV NS3 Biotin labeled protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

HCV NS3 Biotin although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.