

HCV NS4 Mosaic Genotype-1 Recombinant Protein



RPPB5455

Product Information Protein Information

Product SKU:

RPPB5455

Protein description:

The E.coli derived recombinant protein contains the HCV NS4 immunodominant regions, amino acids 1691-1710, 1712-1733 & 1921-1940, fused with a GST tag at N-terminus.

Formulation:

1.5M urea, 25mM Tris-HCl, pH 8.0, 0.2% Triton-X & 50% glycerol.

Purity:

Protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

HCV NS4 Mosaic Genotype-1 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.