

# HCV NS4 Mosaic Genotype-3 Recombinant Protein



RPPB5457

---

## Product Information    Protein Information

**Product SKU:**

RPPB5457

**Protein description:**

The E.coli derived recombinant protein contains the HCV NS4 immunodominant regions, amino acids 1691-1710, 1712-1733 & 1921-1940, fused with a GST tag at N-terminus.

**Formulation:**

1.5M urea, 25mM Tris-HCl, pH 8.0, 0.2% Triton-X & 50% glycerol.

**Purity:**

Protein is >95% pure as determined by 10% PAGE (coomassie staining).

**Stability:**

HCV NS4 Mosaic Genotype-3 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.