

HCV NS4 (1916-1947 a.a.) Recombinant Protein



RPPB5459

Product Information Protein Information

Product SKU:

RPPB5459

Protein description:

E.coli derived 30 kDa recombinant protein. Artificial mosaic polypeptide composite constructed from diagnostically relevant antigenic regions derived from the NS4 region.

Formulation:

1.5M urea, 25mM Tris-HCl pH-8, 0.2% Triton-X & 50% Glycerol.

Purity:

HCV NS4 protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

HCV NS4 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.