

# HCV NS4 Recombinant Protein



RPPB5460

---

## Product Information    Protein Information

**Product SKU:**

RPPB5460

**Host:**

Escherichia Coli.

**Protein description:**

The E.coli derived recombinant Horseradish Peroxidase Labeled protein contains the HCV NS4 immunodominant regions.

**Formulation:**

25mM Tris-HCl pH 8, 1mM EDTA, 1.5M urea & 50% glycerol.

**Purity:**

HCV NS4 Protein is >95% pure as determined by 10% PAGE (coomassie staining).

**Stability:**

HCV NS4 a+b HRP although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.