

HCV NS4 a+b Recombinant Protein



RPPB5464

Product Information Protein Information

Product SKU:

RPPB5464

Protein description:

The E.coli derived 19 kDa recombinant protein Biotin labeled contains the HCV NS4 Genotype 1b immunodominant regions, amino acids 1658-1863. The protein is fused with b-galactosidase (114 kDa) at N-terminus.

Formulation:

(1mg/ml) 20mM Tris-HCl pH 8 and 8M urea.

Purity:

HCV NS4 a+b Biotin protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

HCV NS4 a+b Biotin although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.