

HCV NS5 Genotype-1 Recombinant Protein



RPPB5465

Product Information Protein Information

Product SKU:

RPPB5465

Host:

Escherichia Coli.

Protein description:

The E.coli derived recombinant protein fused to a GST tag contains the HCV NS5 Genotype1 immunodominant regions.

Formulation:

50mM Tris, pH-8 & 5mM EDTA.

Purity:

HCV NS5 Genotype-1 protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

HCV NS5 Genotype-1 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.