

HCV NS5 Genotype 1a Recombinant Protein



RPPB5466

Product Information Protein Information

Product SKU:

RPPB5466

Protein description:

The E.coli derived recombinant protein contains the HCV NS5 Genotype1a immunodominant regions, amino acids 2322-2423 having a total Mw of 38.34kD which includes a 26Kda GST tag.

Formulation:

50mM Tris, pH 8 & 5mM EDTA.

Purity:

HCV NS5 Genotype-1a protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

HCV NS5 Genotype-1a although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.